

INTERFERin[®] transfection reagent

Short protocol – siRNA transfection



Day 0: Cell seeding

→ Seed cells in **V** mL of serum containing medium according to the table below


Quantities per well, dish or flask

Culture vessel	Number of cells*	V = volume of medium during transfection
96-well	2 500 – 7 500	0.2 mL
24-well	15 000 – 35 000	1 mL
12-well	30 000 – 70 000	2 mL
6-well / 35 mm	100 000 – 200 000	4 mL
100 mm / flask 75 cm ²	750 000 – 1.25 x 10 ⁶	15 mL

*For suspension cells, please refer to the complete protocol.

Day 1: Transfection = 1 nM siRNA

→ Transfect cells at **30-50% confluency**



Quantities per well, dish or flask

Culture vessel	W = volume of medium without serum	X = amount of siRNA added (1 nM*)	Y = volume of INTERFERin [®] reagent	Z = volume of growth medium
96-well	50 µL	0.17 pmoles (2.4 ng)	0.75 ± 0.5 µL	0.125 mL
24-well	100 µL	0.6 pmoles (8.4 ng)	2 ± 1 µL	0.5 mL
12-well	200 µL	1.2 pmoles (17 ng)	4 ± 2 µL	1 mL
6-well / 35 mm	200 µL	2.2 pmoles (31 ng)	8 ± 4 µL	2 mL
100 mm / flask 75 cm ²	500 µL	10.5 pmoles (147 ng)	40 ± 10 µL	10 mL

*in final volume of culture.

Day 2-3: Analyze gene silencing

See back page for optimization tips

Download complete protocol on <https://myaccount.polyplus-transfection.com/>

+33 (0)3 90 40 61 80

support@polyplus-transfection.com

www.polyplus-transfection.com

Polyplus

INTERFERin[®] transfection reagent

Short protocol – Optimization tips



+ Protocol Optimization

- + The siRNA final concentration may range from 1 to 50 nM depending on the cells and the target gene.
- + Check our online Cell Transfection Database at:
<http://www.polyplus-transfection.com/resources/cell-transfection-database/>

+ Tips to increase cell viability of sensitive cells

- + Replace medium 4 h after transfection.
- + Check that silencing the target gene does not affect cell viability.

+ Use appropriate controls

- + Positive control: siRNA against housekeeping genes/fluorescently labelled siRNA.
- + Negative control: mismatch, scramble or non-targeting sequence.
- + Be cautious with fluorescently labeled siRNA: 20 to 30 nM are needed to detect a signal, while only 1 nM can be sufficient for efficient silencing using INTERFERin[®].

+ Good siRNA Transfection Practices

- + Store appropriately INTERFERin[®] ($5 \pm 3^\circ\text{C}$). Do NOT freeze INTERFERin[®].
- + Ensure that cells have been passaged more than twice and less than 20 times prior to transfection.
- + Discard overconfluent cells.
- + Check the half-lives of the protein and mRNA of interest and measure gene silencing accordingly 24 to 96 h after transfection.
- + Regularly check for mycoplasma contaminations.
- + Serum quality may drastically affect transfection efficiency. When purchasing a new batch of serum or trypsin, check cell viability as well as transfection efficiency.

Note: INTERFERin[®] is recommended for siRNA transfection. Please refer to the complete protocol available when creating your account online at: <https://myaccount.polyplus-transfection.com/>.

Use jetPRIME[®] for DNA/siRNA co-transfection.